

Installation Instructions for the New Pixera Driver

Camera System: Penguin & Pro Series

Software Applications: Penguin/Pro Application Suite

Operating Systems: MacOS 9.0–9.2

Date: Feb 08 2002

Revision: 1.0

Please, follow the steps below.

Installation of Pixera Driver

Step 1.

Uncompress the driver file, with Stuffit Expander or WinZIP.

Step 2.

Drag & drop the PenguinDriver file into the Extensions folder inside the System Folder on your hard drive. If prompted to replace the current driver, click OK.

Pixera doesn't guarantee compatibility with all vendors but make great effort in providing compatible solutions for a large number of computers, operating systems and other products. In no event will Pixera be Liable for direct, indirect, special, incidental or consequential damages arising out of the use or inability to use the software or documentation, even if advised of the possibility of such damages.

If a recommended solution doesn't work or your system do not comply, please contact our Technical Support at (408) 341 1800 Ext. 3 or via email at techsupport@pixera.com

Copyright © 2002 Pixera Corporation
All rights reserved